

BenzNuclease DNA and RNA Nuclease

Cat. # **BEE-N3116**

For Research Use Only

Description:

BenzNuclease is a recombinant form of *Serratia marcescens* extracellular endonuclease produced in *Escherichia coli* cells using a proprietary process at ACROBiosystems. BenzNuclease is a homodimer with monomer molecular masses about 30 kDa. Two disulfide bonds found in the nuclease are crucial to its activity and stability. The enzyme is a non-specific nuclease with high specific activity, which degrades both single- and double-stranded nucleic acids in any form (single stranded, double stranded, linear, circular and supercoiled). It hydrolyzes internal phosphodiester bonds present between the nucleotides to 5'-phosphorylated oligonucleotides of 3-8 bases in length.

Application :

Its high intrinsic activity and broad substrate tolerance make the endonuclease an ideal tool in a variety of biotechnological and pharmaceutical applications: removal of nucleic acid from protein samples (Elimination of nucleic acids from recombinant proteins; Purification of protein fragments from inclusion bodies; Sample preparation in western blotting or two-dimensional gel electrophoresis) ; Viscosity reduction in protein extracts.

Effect of Conditions on Activity:

BenzNuclease is functional between pH 6 and 10 (optimal at pH8 - 8.5), and from 0°C to 42 °C (optimal at 35 °C - 42 °C). Mg²⁺ (1-2 mM) is required for enzyme activity.

1 mM EDTA reduced the activity by 30% in the presence of 1 mM MgCl₂; 0.1 M EDTA eliminated all enzyme activity. In the presence of 1 mM MgCl₂, enzyme levels were reduced 75% by 0.1 M CaCl₂ or 1 M NaCl. Under standard assay conditions, 1 mM iodoacetate had no effect on the enzymatic rate, whereas 1 mM mercaptoethanol and maleic acid reduced the activity by only 5 to 10%. 10 mM p-Chloromercuribenzoate completely inactivates the enzyme, while 0.64 M beta-mercaptoethanol in the presence of 2 M urea causes only partial inactivation of the enzyme. 4 or 7 M Urea increases the enzyme activity.

Removal of BenzNuclease:

BenzNuclease contain no "Tag" and used in downstream processing can be removed by various purification methods according to the purification strategy for the target protein.

Unit Definition:

One unit will digest sonicated salmon sperm DNA to acid-soluble oligonucleotides equivalent to a ΔA260 of 1.0 in 30 min at pH 8.0 at 37 °C, which corresponds approximately to complete digestion of **37 μg** DNA. Note that 1 KU=1000 units.

Concentration:

≥250 U/μl

Purity:

>99% as determined by SDS-PAGE of reduced BenzNuclease .

SDS-PAGE:




Fig.1 The purity of BenzNuclease was determined by reduced SDS-PAGE and staining overnight with Coomassie Blue.

Formulation:

Lyophilized in Tris HCl, pH 8.0, MgCl₂, and NaCl.

Reconstitution:

See Certificate of Analysis for reconstitution instructions and specific concentrations.

Storage:

Avoid repeated freeze-thaw cycles.

No activity loss was observed after storage at:

In lyophilized state for 1 year (4°C); After reconstitution under sterile conditions for 3 months (-70°C).

Activity Assay Procedure:

1. Reagents and solutions preparation

Reaction buffer*:

50 mM Tris-HCl, 1 mM MgCl₂, pH 8.0 (* In the case of extensive dilution before use, carrier protein such as 0.1 mg/ml HSA or BSA is generally recommended to avoid any enzyme loss from surface adsorption)

DNA Substrate:

1 mg/ml salmon sperm DNA is dissolved overnight at 4 °C, in reaction buffer, and is then sonicated on ice to obtain a homogenous solution.

Enzyme:

Different dilution of nuclease with reaction buffer

Stop reagent:

Trichloroacetic acid (TCA)

2. Standard curve establishment

400 µl substrate + 100 µl enzyme of known activity

= 500 µl mixture

- Incubate the mixture at 37°C for 30 min.
- Stop the reaction by addition of 400 µl cold TCA and incubate on ice for 10 min.
- Centrifuge at 8500 g for 5 min.
- Measure the absorbance of supernatant at 260 nm.
- Lot a standard curve with nuclease of known activities for each set of measurements.

3. Measurement of activity

The activity of any unknown nuclease can be determined from a single measurement by means of the standard curve. The specific activity of BenzNuclease is 1.3×10^6 ~ 1.5×10^6 unit/mg protein.




Fig 2. Standard curve of nuclease activity for BenzNuclease

References:

- (1) Molloy, M.P., et al., 1998, Electrophoresis 19, 837-844.
- (2) Herbert, B.R., 1998, Electrophoresis 19, 845-851.
- (3) Benzonase® Brochure, 1999, Code No. W 220911, Merck KGaA, Darmstadt, Germany.
- (4) Meiss, G., et al., 1995, Biochemistry, 34, 11979–11988

Please contact us via Techsupport@acrobiosystems.com, if you have any question on this product.