



Synonym

ANGPTL7,CDT6

Source

Human ANGPTL7, His Tag(AN7-H52H5) is expressed from human 293 cells (HEK293). It contains AA Gln 27 - Pro 346 (Accession # [O43827-1](#)).

Predicted N-terminus: Gln 27

Molecular Characterization

ANGPTL7(Gln 27 - Pro 346)
O43827-1

Poly-his

This protein carries a polyhistidine tag at the C-terminus.

The protein has a calculated MW of 39.0 kDa. The protein migrates as 32-35 kDa under reducing (R) condition (SDS-PAGE).

Endotoxin

Less than 1.0 EU per µg by the LAL method.

Purity

>85% as determined by SDS-PAGE.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage


For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

SDS-PAGE




Human ANGPTL7, His Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 85%.

Bioactivity-SPR

Discounts, Gifts,
and more!





Human ANGPTL7, His Tag (Cat. No. AN7-H52H5) immobilized on CM5 Chip can bind Human LILRB4, His Tag (Cat. No. LI4-H52H7) with an affinity constant of 4.67 μ M as determined in a SPR assay (Biacore T200) (Routinely tested).

Background

Angiopoietin-like 7 (ANGPTL7) is also known as Corneal-Derived Transcript 6 (CDT6), and is a secreted glycoprotein that is structurally related to the angiopoietins. It is expressed in the corneal stroma, trabecular meshwork, and sclera and is elevated in glaucoma aqueous humor. Recently, ANGPTL proteins, including ANGPTL-2, -3, -4, -5 and -7, are identified as growth factors in combination with SCF, Thrombopoietin, IGF-II and FGF acidic that enhance the expansion and engraftment of human and mouse hematopoietic stem cells.

Clinical and Translational Updates

Please contact us via TechSupport@acrobiosystems.com if you have any question on this product.

Discounts, Gifts,
and more!

