



**Synonym**

CD4,CD4mut,LEU3

**Source**

PE-Labeled Human CD4 Protein, His Tag (CD4-HP2H8) is produced via conjugation of PE to Human CD4 Protein, His Tag with a new generation site-specific technology under Star Staining labeling platform. Human CD4 Protein, His Tag is expressed from human 293 cells (HEK293). It contains AA Lys 26 - Pro 396 (Accession # [AAH25782](#)).

Predicted N-terminus: Lys 26

**Molecular Characterization**

CD4(Lys 26 - Pro 396)  
AAH25782

Poly-his

This protein carries a polyhistidine tag at the C-terminus.

The protein has a calculated MW of 55.8 kDa.

**Conjugate**

PE

Excitation Wavelength: 488 nm / 561 nm

Emission Wavelength: 575 nm

**Purity**

>90% as determined by SDS-PAGE.

**Formulation**

Lyophilized from 0.22 µm filtered solution in PBS, 0.2% BSA, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

**Reconstitution**

Please see Certificate of Analysis for specific instructions.

*For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.*

**Storage**

For long term storage, the product should be stored at lyophilized state at -20°C or lower.

*Please protect from light and avoid repeated freeze-thaw cycles.*

This product is stable after storage at:


- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

**Star Staining** fluorescent-labeled products are developed by a new-generation site-specific labeling technology with Star Standard quality at ACROBiosystems

- ★ Using new-generation site-specific labeling technology to maintain natural bioactivity.
- ★ High specificity and sensitivity verified by flow cytometry.
- ★ No non-specific binding to non-transduced PBMCs.
- ★ High homogeneity and high batch-to-batch consistency.

**Evaluation of CAR expression**

FACS Analysis of Anti-CD4 CAR Expression



Discounts, Gifts, and more!



# PE-Labeled Human CD4 Protein, His TagStar Staining


Catalog # CD4-HP2H8



BIOSYSTEMS  
**Acro**

5e5 of anti-CD4 (Ibalizumab) CAR-293 cells were stained with 100  $\mu$ L of 1:50 dilution (2  $\mu$ L stock solution in 100  $\mu$ L FACS buffer) of PE-Labeled Human CD4 Protein, His Tag (Cat. No. CD4-HP2H8) and negative control protein respectively (Fig. C and B), and non-transfected 293 cells were used as a control (Fig. A). PE signal was used to evaluate the binding activity (QC tested).

FACS Analysis of Non-specific binding to PBMCs



5e5 of PBMCs were stained with PE-Labeled Human CD4 Protein, His Tag (Cat. No. CD4-HP2H8) and anti-CD3 antibody, washed and then analyzed with FACS. FITC signal was used to evaluate the expression of CD3+ T cells in PBMCs, and PE signal was used to evaluate the non-specific binding activity to PBMCs (QC tested).

## Background

CD4 is a glycoprotein that serves as an essential co-receptor on the surface of T lymphocytes (T cells), particularly helper T cells. It plays a critical role in the immune system by recognizing and binding to major histocompatibility complex (MHC) class II molecules on antigen-presenting cells (APCs), thereby facilitating T cell activation and immune response.

Discounts, Gifts,  
and more!

