



Synonym

CEACAM-5,CD66e,CEA,Meconium antigen 100

Source

Cynomolgus CEACAM-5, His Tag(CE5-C52H5) is expressed from human 293 cells (HEK293). It contains AA Gln 35 - Gly 685 (Accession # [XP_005589491.1](#)).

Predicted N-terminus: Gln 35

Molecular Characterization

CEACAM-5(Gln 35 - Gly 685)
XP_005589491.1 Poly-his

This protein carries a polyhistidine tag at the C-terminus.

The protein has a calculated MW of 73.6 kDa. The protein migrates as 100-170 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than 1.0 EU per µg by the LAL method.

Purity

>90% as determined by SDS-PAGE.

>90% as determined by SEC-MALS.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage


For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:


- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

SDS-PAGE



Cynomolgus CEACAM-5, His Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With [Star Ribbon Pre-stained Protein Marker](#)).

SEC-MALS



The purity of Cynomolgus CEACAM-5, His Tag (Cat. No. CE5-C52H5) is more than 90% and the molecular weight of this protein is around 230-260 kDa verified by SEC-MALS.

[Report](#)

Background

Carcinoembryonic antigen-related cell adhesion molecule 5 (CEACAM5) is also Carcinoembryonic antigen (CEA), Meconium antigen 100, CD antigen CD66e, CEACAM5 belongs to the immunoglobulin superfamily and CEA family. CEACAM5 contains seven Ig-like (immunoglobulin-like) domains. CEACAM5 is

Discounts, Gifts,
and more!



Cynomolgus CEACAM-5 / CD66e Protein, His Tag (MALS verified)

Catalog # CE5-C52H5



BIOSYSTEMS
Acro

homodimer protein which binding of E.coli Dr adhesins leads to dissociation of the homodimer. CEACAM5 is cell surface glycoprotein that plays a role in cell adhesion and in intracellular signaling. CEACAM5 is receptor for E.coli Dr adhesins.

Clinical and Translational Updates

Please contact us via TechSupport@acrobiosystems.com if you have any question on this product.

Discounts, Gifts,
and more!



➤ www.acrobiosystems.com

2/4/2024