



Source

Monoclonal Anti-HA (Influenza A/guinea fowl/Hong Kong/WF10/99 (H9N2)) Antibody, Mouse IgG1 (9C1) is a Mouse monoclonal antibody produced from a hybridoma created by fusing SP2/0 myeloma and Mouse B-lymphocytes.

Clone

9C1

Species

Mouse

Isotype

Mouse IgG1 | Mouse Kappa

Conjugate

Unconjugated

Antibody Type

Hybridoma Monoclonal

Reactivity

Virus

Immunogen

Recombinant Influenza A [A/guinea fowl/Hong Kong/WF10/99(H9N2)] Hemagglutinin (HA) Protein is expressed from human 293 cells.

Specificity

Specifically recognizes Influenza A (H9N2) Viruses Hemagglutinin (HA).

Application

Application	Recommended Usage
ELISA	0.2-100 ng/mL

Cross Verification

This product No cross-reactivity in ELISA with
Influenza A [Victoria/4897/2022] Hemagglutinin (HA) Protein, His Tag (Cat. No. HA1-V52H8).
Influenza A [Wisconsin/67/2022] Hemagglutinin (HA) Protein, His Tag (Cat. No. HA1-V52H7).
Influenza A [A/Hong Kong/483/97 (H5N1)] HA, His Tag (Cat. No. HA1-V5229).
Influenza A (Vietnam/1194/2004(H5N1)) Hemagglutinin (HA) Protein, His Tag (Cat. No. HA1-V52H9).
Influenza A (Guangdong/18SF020(H5N6)) Hemagglutinin (HA) Protein, His Tag (Cat. No. HA6-V52H3).
Influenza A (turkey/Germany-MV/R2472/2014(H5N8)) HA Protein, His Tag (Cat. No. HA8-V52H3).
Influenza A (A/Shanghai/02/2013(H7N9)) Hemagglutinin (HA) Protein, His Tag (Cat. No. HA9-V52H3).
This product No cross-reactivity in ELISA with

Purity

>95% as determined by SDS-PAGE.

Purification

Protein A purified/ Protein G purified

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.0 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

Discounts, Gifts,
and more!



Monoclonal Anti-HA (Influenza A/guinea fowl/Hong Kong/WF10/99 (H9N2)) Antibody, Mouse IgG1 (9C1)




Catalog # HA2-Y198

Influenza A Virus (A/District of Columbia/27/2023) HA (H3N2) Protein, His Tag (Cat. No. H32-V52H5).

Influenza A Virus (A/Croatia/10136RV/2023) HA (H3N2) Protein, His Tag (Cat. No. H32-V52H4).


SDS-PAGE



Monoclonal Anti-HA (Influenza A/guinea fowl/Hong Kong/WF10/99 (H9N2)) Antibody, Mouse IgG1 (9C1) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95% (With [Star Ribbon Pre-stained Protein Marker](#)).

Bioactivity-ELISA


Monoclonal Anti-HA (Influenza A/guinea fowl/Hong Kong/WF10/99 (H9N2)) Antibody, Mouse IgG1 (9C1) ELISA
0.2 µg of Influenza A [A/guinea fowl/Hong Kong/WF10/99(H9N2)] Hemagglutinin (HA) Protein, His Tag (MALS verified) per well



Monoclonal Anti-HA (Influenza A/guinea fowl/Hong Kong/WF10/99 (H9N2)) Antibody, Mouse IgG1 (9C1) Conc. (ng/ml)

Immobilized Influenza A [A/guinea fowl/Hong Kong/WF10/99(H9N2)] Hemagglutinin (HA) Protein, His Tag (MALS verified)(Cat. No. HA2-V52H7) at 2 µg/mL (100 µL/well) can bind Monoclonal Anti-HA (Influenza A/guinea fowl/Hong Kong/WF10/99 (H9N2)) Antibody, Mouse IgG1 (9C1)(Cat. No. HA2-Y198) with a linear range of 0.195-6.25 ng/mL (QC tested).

Detection of Monoclonal Anti-HA (Influenza A/guinea fowl/Hong Kong/WF10/99 (H9N2)) Antibody, Mouse IgG1 (9C1) by ELISA Assay



Monoclonal Anti-HA (Influenza A/guinea fowl/Hong Kong/WF10/99 (H9N2)) Antibody, Mouse IgG1 (9C1) Conc. (ng/ml)

Immobilized Influenza A [A/guinea fowl/Hong Kong/WF10/99(H9N2)] Hemagglutinin (HA) Protein, His Tag (MALS verified)(Cat. No. HA2-V52H7) can bind Monoclonal Anti-HA (Influenza A/guinea fowl/Hong Kong/WF10/99 (H9N2)) Antibody, Mouse IgG1 (9C1)(Cat. No. HA2-Y198). The antibody does not bind Influenza A [A/Bangkok/1/1979 (H3N2)] Hemagglutinin (HA) Protein, His Tag (MALS verified)(Cat. No. HA2-V52H3), Influenza A [A/Wisconsin/588/2019 (H1N1)] Hemagglutinin (HA) Protein, His Tag (MALS verified)(Cat. No. HA1-V52H3), Influenza A [A/Darwin/6/2021 (H3N2)] Hemagglutinin (HA) Protein, His Tag (MALS verified)(Cat. No. HA2-V52H5), Influenza A [Sydney/5/2021 (H1N1)] Hemagglutinin (HA) Protein, His Tag (MALS verified)(Cat. No. HA1-V52H4), Influenza B [Austria/1359417/2021 (B/Victoria lineage)] Hemagglutinin (HA) Protein, His Tag (Cat. No. HAE-V52H3), Influenza B [Phuket/3073/2013 (B/Yamagata lineage)] Hemagglutinin (HA) Protein, His Tag (Cat. No. HAE-V52H4), Influenza A [A/Darwin/9/2021 (H3N2)] Hemagglutinin (HA) Protein, His Tag (MALS verified)(Cat. No. HA2-V52H6) and Influenza A [A/Victoria/2570/2019] Hemagglutinin (HA) Protein, His Tag (MALS verified) (Cat. No. HA1-V52H6) (Routinely tested).

Discounts, Gifts,
and more!





Background

Neuraminidase (NA) and hemagglutinin (HA) are major membrane glycoproteins found on the surface of influenza virus. Hemagglutinin binds to the sialic acid-containing receptors on the surface of host cells during initial infection and at the end of an infectious cycle. Hemagglutinin also plays a major role in the determination of host range restriction and virulence. As a class I viral fusion protein, hemagglutinin is responsible for penetration of the virus into the cell cytoplasm by mediating the fusion of the membrane of the endocytosed virus particle with the endosomal membrane.

Clinical and Translational Updates

Discounts, Gifts,
and more!

